

45 min

Hybridization and Bead Preparation

Hybridization

- Prepare 1 - 1,000 ng total RNA in 26 μ l.
- Add 4 μ l **HS** ●.
- Add 5 μ l **PM** ●, mix until homogeneous.
- Denature for 5 min at 75 °C / 1,250 rpm.
- Incubate for 30 min at 60 °C / 1,250 rpm.

Bead Washing

- Resuspend **DB** ●, transfer 75 μ l to a fresh tube.
- Place on magnet for 2 - 5 min, discard supernatant.
- Resuspend beads in 75 μ l **DS** ●, incubate 2 min on magnet, discard supernatant. Repeat once.
- Resuspend beads in 30 μ l **DS** ●.

45 min

Depletion and Purification

Depletion

- Spin down hybridized sample. Add 30 μ l of prepared beads. Mix by pipetting 8x, or until homogeneous.
- Incubate for 15 min at 60 °C / 1,250 rpm. Spin down.
- Place on magnet for 5 min.
- Transfer 60 μ l supernatant to a fresh tube. **ATTENTION:** The supernatant contains the rRNA depleted RNA.

Purification

- Add 24 μ l **PB** ○ and 108 μ l **PS** ○, mix well, incubate for 5 min at RT.
- Place on magnet for 5 - 10 min, discard supernatant.
- Wash the beads twice with 120 - 150 μ l 80 % EtOH, 30 sec. **ATTENTION:** Use 150 μ l for 1.5 ml tubes.
- Air dry beads for 5 - 10 min. **ATTENTION:** do not over dry the beads!
- Add 12 μ l **EB** ○, remove from magnet, mix well, incubate 2 min at RT.
- Place on magnet for 2 - 5 min, transfer 10 μ l of the supernatant to a fresh tube. ⚠ Safe stopping point.

ATTENTION: Spin down solutions before opening tubes or plates!

45 min

Hybridization and Bead Conditioning

Hybridization

- Prepare 1 - 1,000 ng total RNA in 26 μ l.
- Add 4 μ l **HS** ●.
- Add 5 μ l **PM** ●, mix until homogeneous.
- Denature for 5 min at 75 °C / 400 rpm.
- Incubate for 30 min at 60 °C / 400 rpm.

Bead Conditioning

- Resuspend **DB** ●, transfer 75 μ l to a fresh tube.
- Place on magnet for 2 - 5 min, discard supernatant.
- Resuspend beads in 75 μ l **CS** ●, incubate 2 min on magnet, discard supernatant. Repeat once.
- Resuspend beads in 75 μ l **DS** ●, place on magnet for 2 - 5 min, discard supernatant. Repeat twice.
- Resuspend beads in 30 μ l **DS** ●.

75 min

Depletion and Purification

Depletion

- Spin down hybridized sample. Add 30 μ l conditioned beads. Mix by pipetting 8x, or until homogeneous.
- Incubate for 15 min at 60 °C / 400 rpm. Spin down.
- Place on magnet for 5 min.
- Transfer 60 μ l supernatant to a fresh tube. **ATTENTION:** The supernatant contains the rRNA depleted RNA.

Purification

- Add 24 μ l **PB** ○ and 108 μ l **PS** ○, mix well, incubate for 20 min at RT.
- Place on magnet for 5 - 10 min, discard supernatant.
- Add 30 μ l **EB** ○, remove from magnet, mix well, incubate 2 min at RT.
- Add 66 μ l **PS** ○, mix well, incubate 5 min at RT.
- Place on magnet for 2 - 5 min, discard supernatant.
- Wash the beads twice with 120 - 150 μ l 80 % EtOH, 30 sec.
- Air dry beads for 5 - 10 min. **ATTENTION:** do not over dry the beads!
- Add 12 μ l **EB** ○, remove from magnet, mix well, incubate 2 min at RT.
- Place on magnet for 2 - 5 min, transfer 10 μ l of the supernatant to a fresh tube.

ATTENTION: Spin down solutions before opening tubes or plates!