

RT² Profiler PCR Array (96-Well Format and 384-Well [4 x 96] Format)

Mouse MAP Kinase Signaling Pathway

Cat. no. 330231 PAMM-061ZA

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format A	Applied Biosystems® models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad® models iCycler®, iQ™ 5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf® Mastercycler® ep realplex models 2, 2s, 4, 4s; Stratagene® models Mx3005P®, Mx3000P®; Takara TP-800
RT ² Profiler PCR Array, Format C	Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)
RT ² Profiler PCR Array, Format D	Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon®, DNA Engine Opticon 2; Stratagene Mx4000®
RT ² Profiler PCR Array, Format E	Applied Biosystems models 7900HT (384-well block), ViiA 7 (384-well block); Bio-Rad CFX384™
RT ² Profiler PCR Array, Format F	Roche® LightCycler® 480 (96-well block)
RT ² Profiler PCR Array, Format G	Roche LightCycler 480 (384-well block)
RT ² Profiler PCR Array, Format H	Fluidigm® BioMark™



Sample & Assay Technologies

Description

The Mouse MAP Kinase Signaling Pathway RT² Profiler PCR Array profiles the expression of 84 genes related to the MAP kinase (MAPK) signaling pathway. Members of the MKKK, MKK, and MAPK families are represented on this array. Transcription factors and genes whose expression is induced by MAP kinase signaling are included as well. Raf regulating proteins and MEKK1 interacting proteins are found on this array along with cell cycle proteins regulated by the Erk1/2 pathway. Genes related to scaffolding and anchoring are also included. Using real-time PCR, you can easily and reliably analyze expression of a focused panel of genes related to the MAPK signaling pathway with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in formats A, C, D, E, F, and G are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products. RT² Profiler PCR Arrays in format H are shipped on dry ice or blue ice packs.

For long term storage, keep plates at –20°C.

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.

Array layout (96-well)

For 384-well 4 x 96 PCR arrays, genes are present in a staggered format. Refer to the *RT² Profiler PCR Array Handbook* for layout.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Araf	Aif2	Ccna1	Ccna2	Ccnb1	Ccnb2	Ccnd1	Ccnd2	Ccnd3	Ccne1	Cdc42	Cdk2
B	Cdk4	Cdk6	Cdkn1a	Cdkn1b	Cdkn1c	Cdkn2a	Cdkn2b	Cdkn2c	Cdkn2d	Chuk	Col1a1	Creb1
C	Crebbp	Dlk1	E2f1	Egfr	Egr1	Elk1	Ets1	Ets2	Fos	Grb2	Hras1	Hspa5
D	Hspb1	Jun	Kcnn1	Kras	Ksr1	Lamtor3	Map2k1	Map2k2	Map2k3	Map2k4	Map2k5	Map2k6
E	Map2k7	Map3k1	Map3k2	Map3k3	Map3k4	Map4k1	Mapk1	Mapk10	Mapk11	Mapk12	Mapk13	Mapk14
F	Mapk3	Mapk6	Mapk7	Mapk8	Mapk8ip1	Mapk8ip2	Mapk8ip3	Mapk9	Mapkapk2	Mapkapk5	Max	Mef2c
G	Mknk1	Mos	Myc	Nf1atc4	Nras	Pak1	Rac1	Raf1	Rb1	Sfn	Smad4	Trp53
H	Actb	B2m	Gapdh	Gusb	Hsp90ab1	MGDC	RTC	RTC	RTC	PPC	PPC	PPC

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Mm.220946	NM_009703	Araf	V-raf murine sarcoma 3611 viral oncogene homolog
A02	Mm.209903	NM_009715	Aif2	Activating transcription factor 2
A03	Mm.4815	NM_007628	Ccna1	Cyclin A1
A04	Mm.4189	NM_009828	Ccna2	Cyclin A2
A05	Mm.260114	NM_172301	Ccnb1	Cyclin B1
A06	Mm.22592	NM_007630	Ccnb2	Cyclin B2
A07	Mm.273049	NM_007631	Ccnd1	Cyclin D1
A08	Mm.333406	NM_009829	Ccnd2	Cyclin D2
A09	Mm.246520	NM_007632	Ccnd3	Cyclin D3
A10	Mm.16110	NM_007633	Ccne1	Cyclin E1
A11	Mm.1022	NM_009861	Cdc42	Cell division cycle 42 homolog (S. cerevisiae)
A12	Mm.111326	NM_016756	Cdk2	Cyclin-dependent kinase 2
B01	Mm.6839	NM_009870	Cdk4	Cyclin-dependent kinase 4
B02	Mm.31672	NM_009873	Cdk6	Cyclin-dependent kinase 6
B03	Mm.195663	NM_007669	Cdkn1a	Cyclin-dependent kinase inhibitor 1A (P21)
B04	Mm.2958	NM_009875	Cdkn1b	Cyclin-dependent kinase inhibitor 1B
B05	Mm.168789	NM_009876	Cdkn1c	Cyclin-dependent kinase inhibitor 1C (P57)
B06	Mm.4733	NM_009877	Cdkn2a	Cyclin-dependent kinase inhibitor 2A
B07	Mm.423094	NM_007670	Cdkn2b	Cyclin-dependent kinase inhibitor 2B (p15, inhibits CDK4)
B08	Mm.1912	NM_007671	Cdkn2c	Cyclin-dependent kinase inhibitor 2C (p18, inhibits CDK4)
B09	Mm.29020	NM_009878	Cdkn2d	Cyclin-dependent kinase inhibitor 2D (p19, inhibits CDK4)
B10	Mm.3996	NM_007700	Chuk	Conserved helix-loop-helix ubiquitous kinase
B11	Mm.277735	NM_007742	Col1a1	Collagen, type I, alpha 1
B12	Mm.453295	NM_133828	Creb1	CAMP responsive element binding protein 1
C01	Mm.132238	NM_001025432	Crebbp	CREB binding protein
C02	Mm.157069	NM_010052	Dlk1	Delta-like 1 homolog (Drosophila)
C03	Mm.18036	NM_007891	E2f1	E2F transcription factor 1
C04	Mm.8534	NM_007912	Egfr	Epidermal growth factor receptor
C05	Mm.181959	NM_007913	Egr1	Early growth response 1
C06	Mm.405823	NM_007922	Elk1	ELK1, member of ETS oncogene family
C07	Mm.292415	NM_011808	Ets1	E26 avian leukemia oncogene 1, 5' domain
C08	Mm.290207	NM_011809	Ets2	E26 avian leukemia oncogene 2, 3' domain
C09	Mm.246513	NM_010234	Fos	FBJ osteosarcoma oncogene
C10	Mm.439649	NM_008163	Grb2	Growth factor receptor bound protein 2
C11	Mm.334313	NM_008284	Hras1	Harvey rat sarcoma virus oncogene 1
C12	Mm.330160	NM_022310	Hspa5	Heat shock protein 5
D01	Mm.13849	NM_013560	Hspb1	Heat shock protein 1
D02	Mm.275071	NM_010591	Jun	Jun oncogene
D03	Mm.32074	NM_032397	Kcnn1	Potassium intermediate/small conductance calcium-activated channel, subfamily N, member 1
D04	Mm.383182	NM_021284	Kras	V-Ki-ras2 Kirsten rat sarcoma viral oncogene homolog
D05	Mm.4745	NM_013571	Ksr1	Kinase suppressor of ras 1
D06	Mm.331392	NM_019920	Lamtor3	Late endosomal/lysosomal adaptor, MAPK and MTOR activator 3
D07	Mm.248907	NM_008927	Map2k1	Mitogen-activated protein kinase kinase 1
D08	Mm.275436	NM_023138	Map2k2	Mitogen-activated protein kinase kinase 2

Position	UniGene	GenBank	Symbol	Description
D09	Mm.18494	NM_008928	Map2k3	Mitogen-activated protein kinase kinase 3
D10	Mm.412922	NM_009157	Map2k4	Mitogen-activated protein kinase kinase 4
D11	Mm.325746	NM_011840	Map2k5	Mitogen-activated protein kinase kinase 5
D12	Mm.14487	NM_011943	Map2k6	Mitogen-activated protein kinase kinase 6
E01	Mm.3906	NM_011944	Map2k7	Mitogen-activated protein kinase kinase 7
E02	Mm.15918	NM_011945	Map3k1	Mitogen-activated protein kinase kinase kinase 1
E03	Mm.211762	NM_011946	Map3k2	Mitogen-activated protein kinase kinase kinase 2
E04	Mm.27041	NM_011947	Map3k3	Mitogen-activated protein kinase kinase kinase 3
E05	Mm.28587	NM_011948	Map3k4	Mitogen-activated protein kinase kinase kinase 4
E06	Mm.148278	NM_008279	Map4k1	Mitogen-activated protein kinase kinase kinase kinase 1
E07	Mm.196581	NM_011949	Mapk1	Mitogen-activated protein kinase 1
E08	Mm.39253	NM_009158	Mapk10	Mitogen-activated protein kinase 10
E09	Mm.91969	NM_011161	Mapk11	Mitogen-activated protein kinase 11
E10	Mm.38343	NM_013871	Mapk12	Mitogen-activated protein kinase 12
E11	Mm.27970	NM_011950	Mapk13	Mitogen-activated protein kinase 13
E12	Mm.311337	NM_011951	Mapk14	Mitogen-activated protein kinase 14
F01	Mm.8385	NM_011952	Mapk3	Mitogen-activated protein kinase 3
F02	Mm.480076	NM_015806	Mapk6	Mitogen-activated protein kinase 6
F03	Mm.38172	NM_011841	Mapk7	Mitogen-activated protein kinase 7
F04	Mm.21495	NM_016700	Mapk8	Mitogen-activated protein kinase 8
F05	Mm.2720	NM_011162	Mapk8ip1	Mitogen-activated protein kinase 8 interacting protein 1
F06	Mm.173337	NM_021921	Mapk8ip2	Mitogen-activated protein kinase 8 interacting protein 2
F07	Mm.43081	NM_013931	Mapk8ip3	Mitogen-activated protein kinase 8 interacting protein 3
F08	Mm.68933	NM_016961	Mapk9	Mitogen-activated protein kinase 9
F09	Mm.221235	NM_008551	Mapkapk2	MAP kinase-activated protein kinase 2
F10	Mm.272206	NM_010765	Mapkapk5	MAP kinase-activated protein kinase 5
F11	Mm.268548	NM_008558	Max	Max protein
F12	Mm.24001	NM_025282	Mef2c	Myocyte enhancer factor 2C
G01	Mm.209327	NM_021461	Mknk1	MAP kinase-interacting serine/threonine kinase 1
G02	Mm.459300	NM_020021	Mos	Maloney sarcoma oncogene
G03	Mm.2444	NM_010849	Myc	Myelocytomatosis oncogene
G04	Mm.27908	NM_023699	Nfatc4	Nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 4
G05	Mm.400954	NM_010937	Nras	Neuroblastoma ras oncogene
G06	Mm.260227	NM_011035	Pak1	P21 protein (Cdc42/Rac)-activated kinase 1
G07	Mm.292510	NM_009007	Rac1	RAS-related C3 botulinum substrate 1
G08	Mm.184163	NM_029780	Raf1	V-raf-leukemia viral oncogene 1
G09	Mm.273862	NM_009029	Rb1	Retinoblastoma 1
G10	Mm.44482	NM_018754	Sfn	Stratifin
G11	Mm.100399	NM_008540	Smad4	MAD homolog 4 (Drosophila)
G12	Mm.222	NM_011640	Trp53	Transformation related protein 53
H01	Mm.328431	NM_007393	Actb	Actin, beta
H02	Mm.163	NM_009735	B2m	Beta-2 microglobulin
H03	Mm.343110	NM_008084	Gapdh	Glyceraldehyde-3-phosphate dehydrogenase
H04	Mm.3317	NM_010368	Gusb	Glucuronidase, beta
H05	Mm.2180	NM_008302	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
H06	N/A	SA_00106	MGDC	Mouse Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with real-time cyclers that do not require a reference dye, including: Bio-Rad models CFX96, CFX384, DNA Engine Opticon 2; Bio-Rad/MJ Research Chromo4; Roche LightCycler 480 (96-well and 384-well); all other cyclers	330500
RT ² SYBR Green ROX [™] qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Applied Biosystems models 5700, 7000, 7300, 7500 [Standard and FAST], 7700, 7900HT 96-well block [Standard and FAST] and 384-well block, StepOnePlus; Eppendorf Mastercycler ep realplex models 2, 2S, 4, 4S; Stratagene models Mx3000P, Mx3005P, Mx4000; Takara TP-800	330520
RT ² SYBR Green Fluor qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Bio-Rad models iCycler, iQ5, MyiQ, MyiQ2	330510

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group); Applied Biosystems®, ViiA™, StepOnePlus™, ROX™ (Applied Biosystems Corporation or its subsidiaries); Bio-Rad®, iCycler®, iQ™, MyiQ™, Chromo4™, CFX96™, DNA Engine Opticon®, CFX384™ (Bio-Rad Laboratories, Inc.); Stratagene®, Mx3005P®, Mx3000P®, Mx4000® (Stratagene); Eppendorf®, Mastercycler® (Eppendorf AG); Roche®, LightCycler® (Roche Group); Fluidigm® BioMark™ (Fluidigm Corporation); SYBR® (Molecular Probes, Inc.).

1066029 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779



Sample & Assay Technologies